

# Shefaa document :

Base URL : <https://shefaa-backend.vercel.app>

## User end point :

### 1. Forget password

- Method : Post
- /api/auth/forgot-password
- Ex:

```
{
  "identity": "graduationproject911@gmail.com"
}
```

User get a verification code on his email or phone number

### 2. verify reset code

- method : Post
- /api/auth/verify-reset-code
- Ex:

```
{
  "identity": "graduationproject911@gmail.com",
  "code": "1421"
}
```

عشان يتعرف ان ده بتاعه (var احفظولي الايميل او الفون في ) User enter the code he got

### 3. Reset password :

- Method : Post
- /api/auth/reset-password
- Ex:

```
{
  "identity": "graduationproject911@gmail.com",
  "code": "1421",
  "newPassword": "gradProj2026"
}
```

User enter a new password

#### 4. Logout :

- Method: Post
- /api/auth/logout
- Ex:

```
• {  
•   "refreshToken":  
    "eyJhbGciOiJIUzI1NiIsInR5cCI6IkpXVCJ9.eyJpZCI6IjY5OGE1MTEzNTYwMmZlNGZhNmIwND  
    BmYSIsImIhdCI6MTc3MDY3NDI4NiwiZXhwIjoxNzc5Mjc5MDg2fQ.j0SDgXWl09ZeYBPbxK-  
    NG5VhTsMNzpCqEnZtsXXlg4c"  
• }
```

send the refresh token for user who will logout

#### 5. Register :

- Method : post
  - /api/auth/register
1. For patient :

```
2. {  
3.   "name": "grad_parm",  
4.   "username": "gradproj11",  
5.   "email": "gRadAtionPrct929@gmail.com",  
6.   "password": "password123",  
7.   "phoneNumber": "0155555513",  
8.   "role": "patient"  
9. }
```

##### 2. Doctor :

```
3. {  
4.   "name": "grad_parm",  
5.   "username": "gradproj11",  
6.   "email": "gRadAtionPrct929@gmail.com",  
7.   "password": "password123",  
8.   "phoneNumber": "0155555513",  
9.   "role": "doctor"  
10. "membershipPdf"=> PDF file  
11. }
```

Content type => multipart / form data

##### 3. for pharmacy :

| Type        | Value                         | Key                      |
|-------------|-------------------------------|--------------------------|
| Text        | Central Pharmacy              | name                     |
| Text        | central_pharm                 | username                 |
| Text        | pharm@test.com                | email                    |
| Text        | 123456                        | password                 |
| Text        | pharmacy                      | role                     |
| <b>Text</b> | 123456789                     | commercialRegisterNumber |
| <b>File</b> | (اختاري ملف الـ PDF من جهازك) | medicalLicence           |

#### 4. for lab :

| Type        | Value                         | Key                      |
|-------------|-------------------------------|--------------------------|
| Text        | El-Shefaa Lab                 | name                     |
| Text        | shefaalab1                    | username                 |
| Text        | lab@test.com                  | email                    |
| Text        | 123456                        | password                 |
| Text        | lab                           | role                     |
| <b>Text</b> | 987654321                     | commercialRegisterNumber |
| <b>File</b> | (اختاري ملف الـ PDF من جهازك) | medicalLicence           |
| Text        | (أو radiology center) lab     | facilityType             |
| Text        | Dr. Ahmed Ali                 | medicalDirectorName      |
| Text        | DIR-100200                    | directorProfessionalId   |


## 6. Login :

- Method : post
- /api/auth/login
- Ex:


```
{
  "identity": "gRadAtionPrct929@gmail.com",
  "password": "password123"
}
```

## 7. Sign in using google :

- لازم يكون ايميل حقيقي علي جوجل مش ايميل فيك عشان تعرفوا تدخلوا
- بايميل حقيقي postman علي register الاول اعملوا
- url => <https://shefaa-backend.vercel.app/api/auth/google> بعددين خدوا ال
- هيفتح : browser وحطوه في ال



كانكم بتسجلوا بجوجل عادي بقي واختاروا ايميلكم الي لسه مدخليه علي بوستان



المفروض يشتغل ان شاء الله get profile end point => وجربوه علي ال access خدوا ال

## Patient end point :

/\* Bearer token < access token > \*/

1. Get profile :
  - Method : get
  - /api/patient/profile
2. Update basic info :
  - Method : put
  - /api/patient/profile/basic-info
  - Ex:

```
{
  "age": 25,
  "gender": "male",
  "height": 180,
  "weight": 85
}
```

Patient can edit his info

3. Get medical history :

- Method : get
- /api/patient/medical-history

4. Update medical info :

- Method : put
- /api/patient/profile/medical-info
- Ex :

```
{
  "bloodType": "A+",
  "allergies": ["Penicillin", "Dust"],
  "chronicConditions": ["None"]
}
```

5. Add medication :

- Method : post
- /api/patient/medications
- Ex :


```
{
  "name": "panadol",
  "dosage": "500mg",
  "form": "Tablet",
  "timesPerDay": 1,
  "schedule": [ "02:51 PM" ],
  "startDate": "2026-04-24",
  "continuous": true
}
```

السيرفر كل دقيقة هيشوف مين البيشنت الي ليه معاد دوا حالا ويعدهم auto notification \*\* في جزء ال  
console بعدين بيعتلهم ان حالا ليكم دوا اسمه كذا بس دي بتظهر في ال

```
PS D:\fourth year\graduation project\shefaa-backend\Backend> npm run dev
[2026-04-24 14:48:08] INFO: MongoDB Connected: ac-lrlfefk-shard-00-02.m2jkxp8.mongodb.net
Admin account already exists.
[Cron Job] checking at: 02:49 PM
Matching patients found: 0
[Cron Job] checking at: 02:50 PM
Matching patients found: 0
[Cron Job] checking at: 02:51 PM
Matching patients found: 1
[Notification created for: panadol]
[Cron Job] checking at: 02:52 PM
Matching patients found: 0
[Cron Job] checking at: 02:53 PM
Matching patients found: 0
[Cron Job] checking at: 02:54 PM
Matching patients found: 0
[Cron Job] checking at: 02:55 PM
Matching patients found: 0
```

6. Get notification :

- Method : get
- /api/patient/notifications



```
Body Cookies Headers (24) Test Results 200 OK 1.59 s
[
  {
    "_id": "69eb5924046f71ab80f5eb11",
    "recipient": "69d923a5f21cd0e9498c705c",
    "title": "Medication Reminder 🚫",
    "message": "It's time for your panadol dose (500mg).",
    "type": "medication",
    "isRead": false,
    "createdAt": "2026-04-24T11:51:00.330Z",
    "__v": 0
  },
  {
    "_id": "69eb58a3e408c7eb8996b8be",
    "recipient": "69d923a5f21cd0e9498c705c",
    "title": "new medication",
    "message": "you added panadol to your medication list successfully ",
    "type": "medication",
    "isRead": false,
    "createdAt": "2026-04-24T11:48:51.692Z",
    "__v": 0
  },
  {
    "_id": "69eab2b177b2a31e068d682a",
    "recipient": "69d923a5f21cd0e9498c705c",
    "title": "new medication",
    "message": "you added zerc to your medication list successfully ",
    "type": "medication",
```

7. Get medication :

- Method : get
- /api/patient/medications/
- In medication part at UI

8. Confirm medication :

- Method : post
- /api/patient/medications/:medID/confirm
- $\text{confirm}$   $\text{عشان كل ما يضغط علي}$   $\text{adherenceRate}$   $\text{يتحسبه ال}$

9. Notification isRead :

- Method : patch
- /api/patient/notifications/:notificationID/read

#### 10. Appointment upcoming :

- Method : get
- /api/patient/appointments/upcoming

## Doctor end point :

\*\* Bearer token < access token > \*\*

#### 1. Get profile :

- Method : get
- /api/doctor/profile

#### 2. Update profile :

- Method : put
- /api/doctor/profile
- Ex:

```
{
  "specialization": "Cardiology",
  "yearsOfExperience": "10",
  "paymentOption": "in_clinic",
  "about": "Heart specialist"
}
```

#### 3. Add medical record :

- Method : post
- /api/doctor/add-medical-record
- Ex:

```
{
  "patientId": "69875519d4c27be0fe6fbe42",
  "diagnosis": "Severe Bacterial Tonsillitis",
  "prescription": [
    {
      "name": "Augmentin 1g",
      "dosage": "1 tablet",
      "frequency": "Every 12 hours",
      "duration": "7 days"
    },
    {
      "name": "Panadol Advance",
      "dosage": "2 tablets",
      "frequency": "When needed",
      "duration": "3 days"
    }
  ]
}
```



```

• ],
• "notes": "Patient should rest for 3 days.",
• "nextVisitDate": "2026-02-15"}

```

## Lab end points :

**\*\*Bearer token <access token>**

### 1. Add test :

- Method : post
- /api/lab/add-test
- Ex :

```

• {
•   "testName": "Complete Blood Count (CBC)",
•   "price": 200,
•   "estimatedTime": "24 hours"
• }


```

### 2. Get lab's tests :

- Method : get
- /api/lab/my-tests

### 3. Upload result :

- Method : post
- /api/lab/upload-result
- Ex:



none ☒ form-data ☐ x-www-form-urlencoded ☐ raw ☐ binary ☐ GraphQL

| Key                                            | Value                                  | Description | ... | Bulk Edit |
|------------------------------------------------|----------------------------------------|-------------|-----|-----------|
| <input checked="" type="checkbox"/> patientId  | Text 699a2c64df6c8a5e5753c203          |             |     |           |
| <input checked="" type="checkbox"/> testName   | Text CBC                               |             |     |           |
| <input checked="" type="checkbox"/> resultFile | File Chapter 02 - Lexical Analysis.pdf |             |     |           |
| Key                                            | Text Value                             | Description |     |           |

Content type => multipart / form data

4. Get patient result :

- Method : get
- /api/lab/patient-results/:patientId

## Pharmacy end point:

**\*\*Bearer token <access token>\*\***

1. Get stock :

- Method : get
- /api/pharmacy/stock

2. Add medicine :

- Method : post
- /api/pharmacy/add-medicine
- Ex :

```
{
  "medicineName": "ميوفين",
  "quantity": 100,
  "price": 100,
  "expirationDate": "2026-12-01",
  "category": "Other",
  "description": "باسط عضلات"
}
```

3. Update medicine :

- Method : put
- /api/pharmacy/update-medicine/medicineId
- Ex :

```
{
  "quantity": 250
}
```

4. Delete medicine :

- Method : delete
- /api/pharmacy/delete-medicine/medicineId

## 5. Search for medicine :

- Method : get
- /api/pharmacy/search?name=XXX

GET https://shefaa-backend.vercel.app/api/pharmacy/search?name=zertic

Params

| Key  | Value  | Description |
|------|--------|-------------|
| name | zertic |             |

Body

200 OK • 1.24 s • 1.47 KB

```
{
  "data": [
    {
      "_id": "699a269adf6c8a5e5753c1c6",
      "pharmacyId": {
        "_id": "698f27fea0f2e56d5faff5ee"
      },
      "medicineName": "zertic",
      "quantity": 248,
      "price": 100,
      "expirationDate": "2026-12-01T00:00:00.000Z",
      "description": "مضاد هستامين وهو مفيد جدا في حالة الحساسية والحساسية والحساسية",
      "category": "Other",
      "createdAt": "2026-02-21T21:41:46.899Z",
    }
  ]
}
```

## 6. dispense medicine:

- method : post
- /api/pharmacy/dispense
- Ex :

```
{
  "stockId": "699a269adf6c8a5e5753c1c6",
  "quantityToDispense": 2
}
```

## Clinic end point:

### 1. Add new clinic :

- Method : post
- /api/clinic/
- \*\*Bearer token < access token > \*\* => doctor access token
- Ex:

```
{
  "name": "Heart Care Clinic",
}
```

```

• "city": "Cairo",
• "address": "123 Nile St, Maadi",
• "location": { "coordinates": [30.0444, 31.2357] },
• "daysOfWeek": [
•   { "day": "Sunday", "open": "10:00", "close": "18:00" },
•   { "day": "Tuesday", "open": "10:00", "close": "18:00" }
• ],
• "dailyCapacity": 20,
• "slotDuration": 30,
• "capacityPerSlot": 1,
• "price": 500
• }

```

## Appointment end points :

### 1. Search doctor :

- Method : get
- /api/doctor/search-doctors?specialization=cardiology&gender=female&city=cairo

>عشان تدوروا علي الدكتور سواء بالتخصص او مكان العيادة او نوعه بس الاول عدلوا في بروفایل الدكتور عشان تدخلوا نوعه ايه <

### 2. Book appointment :

- Method : post
- \*\* Bearer token <access token > \*\* => patient access token
- /api/appointments/
- Ex:

```

• {
•   "clinicId": "69aaf2fc2b44bad810b01f35",
•   "date": "2026-03-14",
•   "startTime": "5:00 PM",
•   "endTime": "5:30 PM",
•   "appointmentType": "video-call",
•   "paymentOption": "prePay",
•   "isFollowUp": false
• }

```

### 3. Get my appointments ( doctor and patient ):

- Method : get
- \*\* Bearer token <access token > \*\* => علي حسب انتوا دلوقت في بروفایل البيشنت والا الدكتور

- /api/appointments/my

#### 4. Cancel appointment :

- Method : delete
- \*\* Bearer token <access token > \*\* => patient
- /api/appointments/:id
- Id => appointment ID

## Reviews end points :

#### 1. Add review :

- \*\* Bearer token <access token > \*\* => patient
- Method : post
- /api/reviews
- ميقدرش يدي الا ريفيو واحد بس لكل دكتور
- Ex:

```
{
  "doctorId": "698c91524687c848d2b629ea",
  "rating": 3,
  "comment": "she is a good doctor."
}
```

#### 2. get doctor's reviews :

- method : get
- /api/reviews/doctor/:id
- Id => for doctor

#### 3. Get my reviews :

- \*\* Bearer token <access token > \*\* => patient
- Method : get
- /api/reviews/me
- البيشنت بيشوف كل الريفيوز الي عملها

#### 4. Update reviews :

- \*\* Bearer token <access token > \*\* => patient
- Method : put
- /api/reviews/:id

- بتاع الريفير id ال
- Ex:

```

• {
•   "rating": 4,
•   "comment": "تم تعديل التقييم، الخدمة كانت جيدة جداً."
• }

```

#### 5. Delete review :

- \*\* Bearer token <access token > \*\* => patient
- Method : delete
- /api/reviews/:id
- بتاع الريفير id ال

## Admin end points :

- First admin login to get access token :

```

• {
•   "identity": "admin@shefaa.com",
•   "password": "admin123456"
• }

```

#### 1. Get states :

- Method : get
- /api/admin/stats
- Brings the number of registered users in each category

#### 2. get users :

- method : get
- /api/admin/users
- Display a list of all system users, excluding the password.

#### 3. Get user by role :

- Method : get
- /api/admin/users/role/:role

#### 4. Get users pending :

- Method : get
- /api/admin/users/pending
- [ doctor , lab , pharmacy ] profiles which need to activate

#### 5. Activate user :

- Method : patch
- /api/admin/users/activate/:id
- userId => in users database

- activate [ doctor , lab , pharmacy ] profile and automatic send email
6. deactivate user profile :
    - Method : patch
    - /api/admin/users/deactivate/:id
    - userId => in users database
  7. delete user account :
    - method : delete
    - /api/admin/users/:id
    - userId => in users database

## lab report ai end point:

1. upload a pdf or image with any extension of a lab report:
  - method : post
  - /api/labReport/analyze

- Response:

```
{ } JSON ▾ ▶ Preview □ Visualize ▾
1  {
2    "success": true,
3    "data": {
4      "patientName": "سعيدة عزب درويش العبد",
5      "findings": [
6        {
7          "testName": "Erythrocyte Sedimentation Rate (ESR) - First Hour",
8          "result": 50,
9          "unit": "mm",
10         "status": "High"
11       },
12       {
13         "testName": "Erythrocyte Sedimentation Rate (ESR) - Second Hour",
14         "result": 115,
15         "unit": "mm",
16         "status": "High"
17       }
18     ],
19     "dangerScore": 6,
20     "summary": "The ESR levels are significantly elevated, indicating systemic inflammation or infection. However, the current results show a slight improvement compared to the previous readings.",
21     "tips": [
22       "Consult an internal medicine specialist for investigation.",
23       "Stay well hydrated throughout the day.",
24       "Maintain a balanced, anti-inflammatory diet."
25     ]
26   }
27 }
```